


----CHANGE STICKERS EVERYDAY----

Patient Education

- This test is for diagnostic purposes only. In case of an emergency, contact 911 or go to the nearest hospital.
- The cardiac technologists do not monitor patients' data in real time.
- If you experience difficulties please contact our office (905-848-0001 EXT 4).



Guidelines for electrode placement & unit care

- Do not open the monitor at any time.
- Change electrode patches once per day.
- The same electrode patch (sticker) is used for all 3 locations. Move the location at each change within the placement range below.
- **MONITOR IS NOT WATER PROOF.** Remove all Electrode patches before showering or swimming.
- Before putting on electrodes, make sure skin is clean & dry. Do not use any lotion or powder; this is in order to prevent poor electrode contact.
- Do not press or touch the centre portion of the electrode patches where the gel is located. Ensure adhesion of the electrode patches by pressing the edges firmly

----CHANGE STICKERS EVERYDAY----

YOUR APPOINTMENT INFORMATION

- IF YOU ARE UNABLE TO KEEP YOUR APPOINTMENT, PLEASE CONTACT US AS SOON AS POSSIBLE.

DROP OFF APPOINTMENT:

FOLLOW UP APPOINTMENT: